

Customer Information Notice: CIN174

trophon[®]2 Version 1.7 Software Release

Nanosonics has a new software release for the trophon2 device, Version 1.7.

The primary enhancements of this software release include:

- **Added historical cycles look-up tool to reprint the associated labels.**
- **Improved AcuTrace workflow when completing a cycle.**
- **Updated device connectivity functions.**
- **Changed the default *Sleep Timer* settings to OFF.**
- **Added encryption to the *Operator* field when exporting device logs.**
- **Improvements to software security and device reliability.**

The trophon2 User Manual has been updated to reflect the software changes. The latest version of trophon2 User Manual can be accessed online via Nanosonics Academy or your regional Nanosonics Customer Centre.

If your trophon2 device is being upgraded from software versions prior to Version 1.6 you should also refer to the relevant Customer Information Notices (CIN) of previous software releases to ensure that you are familiar with all the changes to your trophon device. The CINs are available in the 'Product & Clinical Information' section of your regional Nanosonics Customer Centre.

The trophon[®] family includes trophon[®] EPR and trophon[®]2 devices which share the same core technology of 'sonically activated' hydrogen peroxide.

Nanosonics, trophon, AcuTrace, NanoNebulant and Sonex-HL are trade marks of Nanosonics Ltd.